

MEDICAL RECOMBINANT HUMAN COLLAGEN REPAIR SPRAY

Long Collagen 1068 contains approximately 200 functional domains, enabling precise, error-free repair. Its Type III collagen is delicate and soft, promoting orderly skin growth while minimizing scarring.

This direct collagen supplement achieves one-step collagen remodeling, effectively repairing skin damage and reducing recovery time.

Benefits:

Fine, Comfortable Mist: The pre-filled metered pump head delivers a gentle, even spray without irritating the skin.

Deep Cleansing & Antibacterial Action: Helps purify the skin while inhibiting bacterial growth.

Non-Contact Application: Ensures safe, hygienic administration without direct skin contact.

Package: 1 box*1.1 fl.oz. /30ml
Shelf life: 2 Years



Effects

General Surgery/ Emergency Department/Burns/Dermatology
Laser Plastic Surgery

Storage

Indoors and transported at room temperature. Low temperature freezing and remelting, does not affect the use.

Key Ingredients

Recombinant human type III collagen, disodium hydrogen phosphate, sodium dihydrogen phosphate and a small amount of sodium chloride, glycerin, water for injection and a small amount of preservatives.

D-NUTRIMEC MECHANISM OF ACTION IN REPAIRING SKIN BARRIER

Barrier

Recombinant Human Type III Collagen forms a reticular structure barrier on the surface of damaged skin to isolate, reduce the invasion of external bacteria and dust.

Guide


Forms a reticular structure at the damaged skin, guides fibroblasts to migrate to the damaged part to divide and proliferate, own collagen synthesis and promotes their own collagen synthesis.

Remodeling

The components work together inside and outside to perfectly play the repair function of active collagen.

Fill

Penetrate into the damaged skin, fill the interstitial gaps, stabilize the skin state, and replenish the missing collagen.





COLLAGEN REPAIR SPRAY

Moisturize and care the skin.

Net. 30ml

BONNEHEURE.



Pioneer of Recombinant Humanized Collagens



BONNEHEURE
Effective Skincare Product

1068

TYPE III COLLAGEN

**Fast
Healing**

**No
Rejection**

**Less
Scarring**

After hyaluronic acid,
a new material for skin repair.



TYPE III COLLAGEN:

quickly repair the skin mucous
membrane barrier function,

ISOLATE ALLERGENS,

form a protective film to lock
the water molecules.




5 BENEFITS

15S

Stop the Bleeding

- ▶ Stops bleeding and pain
- ▶ Supports the skin Scar-Free Repair (Type III Collagen)
- ▶ Moisturizing and water locking
Brightens skin
(lightens discoloration)



Barrier ▸ Fill ▸ Lead ▸ Reshape

• Barrier

Recombinant human collagen and hyaluronic acid form a mesh structure on the mucosal surface to lock in water and moisturize.

• Fill


When the skin breaks down or atrophies, the cellular gaps in the stratum corneum increase and collagen is able to penetrate the stratum corneum.

• Lead

Replenish the defective Chuan collagen in the dermis to repair the damage while promoting the division and proliferation of fibroblasts after they migrate to the damaged area.

• Reshape

The internal and external synergies of the components work perfectly to bring out the function of active collagen and achieve collagen remodeling.



Collagen Skin Care

Increased Skin Elasticity After Application



- ✓ Robust skin barrier
- ✓ Repairing damaged skin
- ✓ Increased elasticity, tenderness and smoothness

28-day New Collagen-empowered Skin, Validated By Authoritative Institutions

28-75 days = a skin growth cycle

Overall recognition of skin

state improvement: 93.3%


+23.92%

Rehydration and
water-locking
Cuticle water content



-21.98%

Barrier function
Trans epidermal
water loss value




-13.80%

Redness reduction effect
a value in redness zone



USE EFFECT



SCOPE OF APPLICATION



Laser plastic surgery

- >Repairing photoelectric damage
- >Rapid relieving redness
- >Reducing blackening after surgery
- >Locking moisture



Surgery trauma

- >Hemostasis and relieving pain
- >Rapidly healing wound
- >Reducing scarring
- >Avoiding pigmentation



Skin diseases

- >Restoring skin function
- >Reducing inflammatory reaction
- >Thickening the cuticle
- >Adjusting water oil balance
- >Repair Acne Damage
- >Avoid acne pits and marks



Current maintenance

- >Supplementing collagen
- >Reducing wrinkle and skin laxity
- >Improving skin texture and complexion
- >Locking skin moisture
- >After Sun Relief
- >Soothes skin after shaving/tattooing

**Beauty Repair ALL
Healthy Ingredients
9 items 0 additions**

0 ADDITIVE

- ▶ Flavoring
- ▶ Coloring
- ▶ Alcohol
- ▶ Salicylic Acid
- ▶ Hormones
- ▶ Retinol
- ▶ Mineral Oil
- ▶ Heavy Metals
- ▶ Formaldehyde Release



SERIES PRODUCTS

Project 28-day rate of change	Serum Gel	Repair Spray	Repair Cream	Repair Mask
				
Skin stratum corneum moisture content (Water locking function)	+23.92%	+23.92%	+23.92%	+23.92%
Transepidermal water loss value (Barrier function)	-21.98%	-14.99%	-15.21%	-16.82%
Skin red zone a* value (Redness reduction and sensitization)	-13.80%	-9.90%	-13.09%	-20.22%